

Communication with technical tools



- **Portable devices of communication**
- **Stationary system:** a personal computer, special communicative software, the device input (output) voice message (touch screen, foot-mouse, pointer-helmet, alternative keyboard and etc.).



Auxiliary devices for use of the technical tools of communication



GoTalk Button




Record and play a 10 second message with this handy and versatile device.

You can use this **for all sorts of scheduling activities** to foster communication and encourage daily living skills.

There's a built-in magnet for display on any metal surface.

There's a record lock feature to eliminate accidental erasures.

BIGmack



Record any single message up to two minutes in length directly into the BIGmack communicator and press its activation surface for playback.

Connect a toy or appliance for additional motivation.

With its large, **5-in/12.7-cm activation surface**, BIGmack is an excellent choice for persons with visual impairments and persons with physical disabilities who require a larger target area.

Powered by a replaceable **9-volt battery**, BIGmack never needs charging and can last for weeks or even months on a single battery.

LITTLEmack



Simply record any message directly into the LITTLEmack communicator for one touch, single message playback up to **two minutes in length**.

Connect a toy or battery-operated appliance for instant positive reinforcement.

With its angled base and **2.5-in/6.4-cm** activation surface, LITTLEmack is perfect for mounting and is recommended for persons who can access a smaller target area.

Powered by a replaceable **9-volt battery**, LITTLEmack never needs charging and can last for weeks or even months on a single battery.

Personal Talker



The pocket-sized **Personal Talker** has **ten seconds of recording time**, enough time to record several sentences.

A pen point is needed to record, eliminating accidental erasure/recording.

Place a picture symbol on the front cover, then simply press on the cover to activate the message.

Talkbook Four



Each GoTalk One is a **high-quality single-message talker** with a large picture display and a **big, colorful PLAY button**.

Put your picture display into the covered frame and slide the switch on the back to RECORD.

Press the PLAY button and talk for up to **10 seconds**. Then slide the switch and PLAY your message.

A durable, affordable, easy-to-carry, four-message communication device.

Go Talk Pocket



Six message keys with five levels give the user plenty to talk about.

Overlays slide in easily and are stored in a removable compartment on the back.

The **GoTalk Pocket** is rugged, attractive, easy to use and has great sound quality.

Size: 3½ x 5½ x 1 inches.

Weight: 7 oz. **Total recording time:** 5 minutes.

Includes three AAA batteries and 2-year warranty.

Go Talk 4+



GoTalk 4+ has a **20 message** capacity (4 large keys, each 3" square, and five recording levels). Plus **two core messages** that remain constant when you change levels.

The big (3" square) message keys let you use large picture cues. And the entire message area is active, so users can touch any part of the image to communicate. A light touch is all that's needed.

Place photos or illustrations on the blank overlays (included) to indicate your recorded message.


Total recording time is 4.25 minutes.

Size: 9 x 12 x 1 $\frac{1}{8}$ inches.

Weight: 23 oz.

Includes Two AA batteries, printed **User's Guide** and a 2-year warranty.

Go Talk 9+



GoTalks provide fantastic sound with volume control, easy sequential recording, quick level erase, built-in overlay storage, record lock, level lock and 2-year warranty.


Features: 9 message keys, 5 recording levels, 45 message capacity, 8.25 minutes of recording time.

Size: 9 x 12 x 1 1/8 inches.

Weight: 23 oz.

Two AA batteries and printed User's Guide included.

Go Talk 20+




Go Talk 20+ has a **100 message** capacity.

Each of the 20 message keys has **5 levels** with **5 buttons** for essential messages which stay the same even when the level is changed.

Overlays are stored inside and slide in and out easily.

Requires 2 AA batteries.

Talking Photo Album



The Albums can be applied to activities for **children or adults with special needs** or for anyone who just wishes to have audible information readily available as a quick reference.

20 pages in a 'flip style' Album. Each of the 20 pages will hold a single 5x7" Photo or Image, (125x175mm).

Record a **separate voice message** on each page using the built-in microphone. Total recording time: **200 minutes**, across 20 pages (e.g. approx 10 minutes per page!)

Size: 225 x 175 x 35mm Approx.

Black Faux Leather Cover.

Easy to use.

SuperTalker



- **Four different grid formats** with matching keyguards
- **16 minutes of recording time** with variable message lengths
- Compatible with Matrix Maker or BoardMaker Overlay software
- Built-in keyguard storage
- **8 levels**
- **8 input jacks** for single-switch access to any specific location
- **Step-ahead jack** for single switch sequential message playback
- **2 Toy/Appliance output jacks**
- **Volume control** with lockout feature
- **Crystal clear sound**

QuickTalker Freestyle



The QuickTalker Freestyle is a flexible, dynamic display **speech generating device (SGD)** with the critical options necessary for speech-language pathologists to better meet the unique needs of their clients with communication disabilities.

QuickTalker Freestyle includes:

- 9.7-in dynamic multi-touch display
- Amplified sound output
- Blue2 alternative access switch
- 2-year warranty